

External dimensions

Height (mm)	210
Width (mm)	250
Depth(mm)	150

Internal dimensions

Height (mm)	130
Width (mm)	190
Depth (mm)	90
Storage height (mm)	70
Storage width (mm)	180
Storage depth (mm)	50

Volumes

Gross volume (litres)	2.2
Storage volume (litres)	-

Ship dimensions

Weight net (kg)	1.2
Height (mm) Cardboard	230
Width (mm) Cardboard	285
Depth (mm) Cardboard	160
Volume (m3) Cardboard	0.01
Weight (kg) Cardboard	2
Height (mm) Wood	-
Width (mm) Wood	-
Depth (mm) Wood	-
Volume (m3) Wood	-
Weight (kg) Wood	-

Accessories

2 Ice Packs of 0.3L crystal thermometer insert
--

Option

Testo Logger 174T/ starter set or LogTag TRID30-7R or LogTag TRED30-7R
--

Passive

Vaccine transport box

RCW2

Technical details

Colour	blue according to RAL5015
CFC/HCFC free	yes
External / internal material	Polyethylene
Insulation material	30 mm Polyurethane foam / Foaming agent; Cyclopentane
Weight fully loaded (kg)	2.1
Required Ice Packs	2 x 0.3L
Required Eutetics	-

Performance

Hold over time at -10°C	-
Hold over time at +22°C	-
Hold over time at +32°C	-
Hold over time at 43°C	-
Cold life without lid openings +2°C to +10°C at ambient 32°C	-
Cold life without lid openings +2°C to +10°C at ambient 43°C	-
Cold life without lid openings -3°C to +10°C	22h at 32°C / 17h at 43°C
Cold life at +43°C	-
Warm life at -20°C	-
Cool life at +43°C	-
Vaccine storage capacity (litres)	0.6
Climate Zone	-

Certification

WHO test procedure	E4 / Proc / 1
WHO specifications	E4 / VC.1
WHO test report	Univalle E4/2009
PQS Code	-
PIS Code	E4/53-M
ADR/RID MARITIM IATA / ICAO	UN 4H2 / X 2.3 / S / "YYYY" B / 921-080151