

External dimensions

Height (mm)	437
Width (mm)	588
Depth(mm)	288

Internal dimensions

Height (mm)	245
Width (mm)	460
Depth (mm)	180
Storage height (mm)	202
Storage width (mm)	326
Storage depth (mm)	107

Volumes

Gross volume (litres)	20
Storage volume (litres)	7.05

Ship dimensions

Weight net (kg)	6.8
Height (mm) Cardboard	470
Width (mm) Cardboard	610
Depth (mm) Cardboard	310
Volume (m3) Cardboard	0.089
Weight (kg) Cardboard	11
Height (mm) Wood	-
Width (mm) Wood	-
Depth (mm) Wood	-
Volume (m3) Wood	-
Weight (kg) Wood	-

Accessories

10 Ice Packs of 0.6L 2 Ice Packs of 0.3L crystal thermometer
--

Option

Testo Logger 174T/ starter set or LogTag TRID30-7R or LogTag TRED30-7R
--

Passive Vaccine transport box RCW8

Technical details

Colour	blue according to RAL5015
CFC/HCFC free	yes
External / internal material	Polyethylene
Insulation material	50-60 mm Polyurethane foam / Foaming agent; Cyclopentane
Weight fully loaded (kg)	16.4
Required Ice Packs	10 x 0.6L and 2 x 0.3L
Required Eutetics	-

Performance

Hold over time at -10°C	-
Hold over time at +22°C	-
Hold over time at +32°C	-
Hold over time at 43°C	-
Cold life without lid openings +2°C to +10°C at ambient 32°C	-
Cold life without lid openings +2°C to +10°C at ambient 43°C	-
Cold life without lid openings -3°C to +10°C	-
Cold life at +43°C	57.9h
Warm life at -20°C	21.6h
Cool life at +43°C	12h
Vaccine storage capacity (litres)	6
Climate Zone	-

Certification

WHO test procedure	E4 / CB01-VP.1
WHO specifications	E4 / CB01.1
WHO test report	Intertek A65524b
PQS Code	E004/003
PIS Code	-
ADR/RID MARITIM IATA / ICAO	UN 4H2 / Y19 / S / "YYYY" B / 921-110058