

External dimensions

Height (mm)	499
Width (mm)	550
Depth(mm)	475

Internal dimensions

Height (mm)	270
Width (mm)	340
Depth (mm)	260
Storage height (mm)	209
Storage width (mm)	251
Storage depth (mm)	176

Volumes

Gross volume (litres)	24
Storage volume (litres)	9.23

Ship dimensions

Weight net (kg)	11.7
Height (mm) Cardboard	530
Width (mm) Cardboard	570
Depth (mm) Cardboard	490
Volume (m3) Cardboard	0.148
Weight (kg) Cardboard	17
Height (mm) Wood	-
Width (mm) Wood	-
Depth (mm) Wood	-
Volume (m3) Wood	-
Weight (kg) Wood	-

Accessories

14 Ice Packs of 0.6L crystal thermometer
--

Option

Testo Logger 174T/ starter set or LogTag TRID30-7R or LogTag TRED30-7R
--

Passive Vaccine transport box RCW12

Technical details

Colour	blue according to RAL5015
CFC/HCFC free	yes
External / internal material	Polyethylene
Insulation material	90-105 mm Polyurethane foam / Foaming agent; Cyclopentane
Weight fully loaded (kg)	23.3
Required Ice Packs	14 x 0.6L
Required Eutetics	-

Performance

Hold over time at -10°C	-
Hold over time at +22°C	-
Hold over time at +32°C	-
Hold over time at 43°C	-
Cold life without lid openings +2°C to +10°C at ambient 32°C	-
Cold life without lid openings +2°C to +10°C at ambient 43°C	-
Cold life without lid openings -3°C to +10°C	-
Cold life at +43°C	114.9h
Warm life at -20°C	40.9h
Cool life at +43°C	26.4h
Vaccine storage capacity (litres)	7
Climate Zone	-

Certification

WHO test procedure	E4 / CB01-VP.1
WHO specifications	E4 / CB01.1
WHO test report	Intertek A65524c
PQS Code	E004/004
PIS Code	-
ADR/RID MARITIM IATA / ICAO	UN 4H2 / X 29.6 / S / "YYYY" B / 921-080154